

REQUEST FOR QUOTATION NO. RFQ24/02854

PURCHASE OF CRIME SCENE INVESTIGATION KITS FOR THE POLICE ACADEMY'S FORENSIC LABORATORY

SECTION 1 REQUEST FOR QUOTATION (RFQ)

UNDP through **Strengthening the Police Training Capacities in the Republic of Moldova Project** kindly requests your quotation for the provision of **RFQ24/02854: Purchase of Crime Scene Investigation Kits for the Police Academy's Forensic Laboratory** as detailed in line items section of this RFQ.

This Request for Quotation comprises the following documents.

Section 1 This RFQ document generated by the online system

Section 2 RFQ Instructions and Data

Annex 1 Schedule of Requirements

Annex 2 Quotation Submission Form

Annex 3 Technical and Financial Offer

Annex 4 Technical Responsiveness Table

When preparing your quotation, please be guided by the RFQ Instructions and Data. Please note that quotations must be submitted directly in the system responding to the questions and uploading required documents by the date and time indicated in the online portal. It is your responsibility to ensure that your quotation is submitted before the deadline. Quotations received after the submission deadline outside the online portal, for whatever reason, will not be considered for evaluation.

Quotations must be submitted directly in Quantum NextGenERP supplier portal following this link <http://supplier.quantum.partneragencies.org> using the profile you may have in the portal (please log in using your username and password).

Follow the instructions in the user guide to search for the tender using search filters, namely **Negotiation ID UNDP-MDA-00446** and subscribe to the tender in order to get notifications in case of amendments of the tender document and requirements.

In case you have never registered before, follow this link to register a profile: <https://estm.fa.em2.oraclecloud.com/fscmUI/faces/PrcPosRegisterSupplier?prcBuld=300000127714247>

Please note that the access link to the Supplier registered profile is sent from Oracle within up to 3 days. In case you have not received the access link after 3 days since registration, you should address for support to UNDP at the email address: sc.md@undp.org. In case you encounter errors with registration (e.g. system states Supplier already is registered), you should address for support to UNDP at the email address: sc.md@undp.org.

Computer firewall could block *oracle* or *undp.org extension* and Suppliers might not receive the Oracle notifications. Please turn down any firewalls on your computers to ensure receipt of email notification.

Do not create a new profile if you already have one. Use the forgotten password feature in case you do not remember the password or the username from previous registration.

Should you require further clarifications on the application through the Quantum online portal, kindly contact the Procurement Unit at sc.md@undp.org. Please pay attention that the bid shall be submitted online through the Quantum system and any bid sent to the above email shall be disqualified.

Should you require further clarifications on the Request for Quotation, Terms of Reference or other requirements, kindly communicate using the messaging functionality in the portal.

Deadline for Submission of Offers (Date and Time), which is visible in the online procurement system will be final. System will not accept submission of any bid after that date and time. It is the responsibility of the bidder to make sure that the bid is submitted prior to this deadline for submission.

Bidders are advised to upload bid documents and to submit their offer a day prior or well before the date and time indicated under the deadline for submission of Offers. Do not wait until last minute. If Bidder faces any issue during submitting offers at the last minutes prior to the deadline for submission, UNDP may not be able to assist on such a short notice and will not be held liable in such instance. UNDP will not accept any offer that is not submitted directly through the System.

Thank you and we look forward to receiving your quotation.

UNDP Moldova

SECTION 2 GENERAL INSTRUCTIONS

Introduction	<p>Bidders shall adhere to all the requirements of this RFQ, including any amendments made in writing by UNDP. This RFQ is conducted in accordance with the UNDP Programme and Operations Policies and Procedures (POPP) on Contracts and Procurement</p> <p>Any Bid submitted will be regarded as an offer by the Bidder and does not constitute or imply the acceptance of the Bid by UNDP. UNDP is under no obligation to award a contract to any Bidder as a result of this RFQ.</p> <p>UNDP reserves the right to cancel the procurement process at any stage without any liability of any kind for UNDP, upon notice to the bidders or cancellation of the tender in the online portal.</p>
Deadline for the Submission of Quotation	<p>Deadline is indicated in the online portal.</p> <p>If any doubt exists as to the time zone in which the quotation should be submitted, refer to http://www.timeanddate.com/worldclock/.</p>
Method of Submission	<p>Quotations must be submitted as follows</p> <p>NextGenERP supplier portal following this link http://supplier.nextgenerp.partneragencies.org/ using the profile you may have in the portal.</p> <p>Follow the instructions in the user guide to search for the tender using Negotiation ID. In case you have never registered before, follow this link to register a profile</p> <p>https://estm.fa.em2.oraclecloud.com/fscmUI/faces/PrcPosRegisterSupplier?prcBuld=300000127715297&_adf.ctrl-state=8godmwdd9_239&_afLoop=7321111756612874&_afWindowMode=0&_afWindowId=null&_af_rFS=16&_afMT=screen&_afMFW=1920&_afMFH=880&_afMFDW=1920&_afMFDH=1080&_afMFC=8&_afMFCI=0&_afMFM=0&_afMFR=96&_afMFG=0&_afMFS=0&_afMFO=0</p> <p>Do not create a new profile if you already have one. Use the forgotten password feature in case you do not remember the password or the username from previous registration.</p> <ul style="list-style-type: none"> File Format All attachments must be in PDF format unless otherwise instructed by UNDP. File names must be in Latin alphabet/keyboard and clearly indicate the content of the document to facilitate review. All files must be free of viruses and not corrupted.
Cost of preparation of quotation	<p>UNDP shall not be responsible for any costs associated with a Supplier's preparation and submission of a quotation, regardless of the outcome or the manner of conducting the selection process.</p>
Supplier Code of Conduct, Fraud, Corruption,	<p>All prospective suppliers must read the United Nations Supplier Code of Conduct and acknowledge that it provides the minimum standards expected of suppliers to the UN. The Code of Conduct, which includes principles on labour, human rights, environment and ethical conduct may be found at https://www.un.org/Depts/ptd/about-us/un-supplier-code-conduct</p> <p>Moreover, UNDP strictly enforces a policy of zero tolerance on proscribed practices, including fraud, corruption, collusion, unethical or unprofessional practices, and obstruction of UNDP vendors and requires all bidders/vendors to observe the highest standard of ethics during the procurement process and contract implementation. UNDP's Anti-Fraud Policy can be found at http://www.undp.org/content/undp/en/home/operations/accountability/audit/office_of_audit_and_investigation.html#anti</p>

Gifts and Hospitality	Bidders/vendors shall not offer gifts or hospitality of any kind to UNDP staff members including recreational trips to sporting or cultural events, theme parks or offers of holidays, transportation, or invitations to extravagant lunches, dinners or similar. In pursuance of this policy, UNDP (a) Shall reject a bid if it determines that the selected bidder has engaged in any corrupt or fraudulent practices in competing for the contract in question; (b) Shall declare a vendor ineligible, either indefinitely or for a stated period, to be awarded a contract if at any time it determines that the vendor has engaged in any corrupt or fraudulent practices in competing for, or in executing a UNDP contract.
Conflict of Interest	<p>UNDP requires every prospective Supplier to avoid and prevent conflicts of interest, by disclosing to UNDP if you, or any of your affiliates or personnel, were involved in the preparation of the requirements, design, specifications, cost estimates, and other information used in this RFQ. Bidders shall strictly avoid conflicts with other assignments or their own interests, and act without consideration for future work. Bidders found to have a conflict of interest shall be disqualified.</p> <p>Bidders must disclose in their Bid their knowledge of the following a) If the owners, part-owners, officers, directors, controlling shareholders, of the bidding entity or key personnel who are family members of UNDP staff involved in the procurement functions and/or the Government of the country or any Implementing Partner receiving goods and/or services under this RFQ.</p> <p>The eligibility of Bidders that are wholly or partly owned by the Government shall be subject to UNDP's further evaluation and review of various factors such as being registered, operated and managed as an independent business entity, the extent of Government ownership/share, receipt of subsidies, mandate and access to information in relation to this RFQ, among others. Conditions that may lead to undue advantage against other Bidders may result in the eventual rejection of the Bid.</p>
Currency of Quotation	<p>Prices shall be quoted only in the currency indicated in the system</p> <p>MDL (Moldovan Leu) for local suppliers</p> <p>and</p> <p>USD (US Dollars) for international suppliers.</p> <p>For evaluation purposes, bids submitted in any other currencies shall be recalculated at UN Operational Rate of Exchange on the submission deadline</p> <p>https://treasury.un.org/operationalrates/OperationalRates.php</p> <p>UNDP shall not be kept liable for any fluctuations of the exchange market during contract implementation, the Contractor being legally responsible to register any loss/gain of currency exchange resulting from payments against the Contract in accordance with the national legislation.</p>
Joint Venture, Consortium or Association	<p>If the Bidder is a group of legal entities that will form or have formed a Joint Venture (JV), Consortium or Association for the Bid, they shall confirm in their Bid that (i) they have designated one party to act as a lead entity, duly vested with authority to legally bind the members of the JV, Consortium or Association jointly and severally, which shall be evidenced by a duly notarized Agreement among the legal entities, and submitted with the Bid; and (ii) if they are awarded the contract, the contract shall be entered into, by and between UNDP and the designated lead entity, who shall be acting for and on behalf of all the member entities comprising the joint venture, Consortium or Association. Refer to Clauses 19 – 24 under Solicitation policy for details on the applicable provisions on Joint Ventures, Consortium or Association.</p>

Only one Bid	<p>The Bidder (including the Lead Entity on behalf of the individual members of any Joint Venture, Consortium or Association) shall submit only one Bid, either in its own name or, if a joint venture, Consortium or Association, as the lead entity of such Joint Venture, Consortium or Association.</p> <p>Bids submitted by two (2) or more Bidders shall all be rejected if they are found to have any of the following</p> <ul style="list-style-type: none"> a) they have at least one controlling partner, director or shareholder in common; or b) any one of them receive or have received any direct or indirect subsidy from the other/s; or b) they have the same legal representative for purposes of this RFQ; or c) they have a relationship with each other, directly or through common third parties, that puts them in a position to have access to information about, or influence on the Bid of, another Bidder regarding this RFQ process; d) they are subcontractors to each other's Bid, or a subcontractor to one Bid also submits another Bid under its name as lead Bidder; or e) some key personnel proposed to be in the team of one Bidder participates in more than one Bid received for this RFQ process. This condition relating to the personnel, does not apply to subcontractors being included in more than one Bid.
Price variation	No price variation due to escalation, inflation, fluctuation in exchange rates, or any other market factors shall be accepted at any time during the validity of the quotation after the quotation has been received.
Alternative Quotes	If alternative quote is permitted, it may be submitted only if a conforming quote to the RFQ requirements is submitted. Where the conditions for its acceptance are met, or justifications are clearly established, UNDP reserves the right to award a contract based on an alternative quote. If multiple/alternative quotes are being submitted, they must be clearly marked as "Main Quote" and "Alternative Quote" directly in the portal and in any supporting document as relevant.
Contact Person for correspondence, notifications and clarifications	<p>Must be submitted directly in the portal using the messaging functionality.</p> <p>Any delay in UNDP's response shall not be used as a reason for extending the deadline for submission, unless UNDP determines that such an extension is necessary and communicates a new deadline to the Proposers.</p>
Right not to accept any quotation	UNDP is not bound to accept any quotation, nor award a contract or Purchase Order
Right to vary requirement at time of award	At the time of award of Contract or Purchase Order, UNDP reserves the right to vary (increase or decrease) the quantity of services and/or goods, by up to a maximum twenty-five per cent (25%) of the total offer, without any change in the unit price or other terms and conditions.
Publication of Contract Award	UNDP will publish the contract awards on the websites of the CO and the corporate UNDP Web site.
Policies and procedures	This RFQ is conducted in accordance with UNDP Programme and Operations Policies and Procedures
UNGM registration	Any Contract resulting from this RFQ exercise will be subject to the supplier being registered at the appropriate level on the United Nations Global Marketplace (UNGM) website at www.ungm.org . The Bidder may still submit a quotation even if not registered with the UNGM, however, if the Bidder is selected for Contract award, the Bidder must register on the UNGM prior to contract signature.

SECTION 2 SPECIAL INSTRUCTIONS

General Conditions of Contract	Any Purchase Order or contract that will be issued as a result of this RFQ shall be subject to one of the General Conditions of Contract below as applicable in each case specified in the Requirements section Applicable GTC <input checked="" type="checkbox"/> General Terms and Conditions / Special Conditions for Contract. Applicable Terms and Conditions and other provisions are available at UNDP/How-we-buy
Special Conditions of Contract	<input checked="" type="checkbox"/> Cancellation of PO/Contract if the delivery/completion is delayed by 30 calendar days <input checked="" type="checkbox"/> Liquidates damages shall be imposed as follows: - 0.3% of contract price per day of delay, up to a maximum of 30 calendar days, after which UNDP may terminate the contract.
Duties and taxes	Article II, Section 7, of the Convention on the Privileges and Immunities provides, inter alia, that the United Nations, including UNDP as a subsidiary organ of the General Assembly of the United Nations, is exempt from all direct taxes, except charges for public utility services, and is exempt from customs restrictions, duties, and charges of a similar nature in respect of articles imported or exported for its official use. All quotations shall be submitted net of any direct taxes and any other taxes and duties, unless otherwise specified in the requirements section. All prices must <input checked="" type="checkbox"/> be exclusive of VAT and other applicable indirect taxes
Eligibility	A vendor who will be engaged by UNDP may not be suspended, debarred, or otherwise identified as ineligible by any UN Organization or the World Bank Group or any other international Organization. Vendors are therefore required to disclose to UNDP whether they are subject to any sanction or temporary suspension imposed by these organizations. Failure to do so may result in termination of any contract or PO subsequently issued to the vendor by UNDP. It is the Bidder's responsibility to ensure that its ultimate beneficial owners, employees, joint venture members, sub-contractors, service providers, suppliers and/or their employees meet the eligibility requirements as established by UNDP. Bidders must have the legal capacity to enter a binding contract with UNDP and to deliver in the country, or through an authorized representative.
Language of quotation	English or Romanian Including supporting documentation as applicable.
Quotation validity period	Quotations shall remain valid for 90 days from the deadline for the Submission of Quotation.
Partial Quotes	<input checked="" type="checkbox"/> Not Permitted
Alternative Quotes	<input checked="" type="checkbox"/> Not permitted
Payment Terms	<input checked="" type="checkbox"/> 100% within 30 days after receipt of goods, works and/or services and submission of payment documentation.
Conditions for Release of Payment	<input checked="" type="checkbox"/> Passing Inspection, Complete Installation (where applicable) <input checked="" type="checkbox"/> Written Acceptance of Goods based on full compliance with RFQ requirements; <input checked="" type="checkbox"/> Complete delivery of contracted goods;
Clarifications	Bidders must send their inquiries and requests for clarifications using the messaging functionality in the portal. <u>PLEASE PAY ATTENTION</u> QUOTES SHALL NOT BE SUBMITTED TO ANY EMAIL ADDRESS BUT ONLY THROUGH THE PORTAL. Requests for clarification from bidders will not be accepted any later than 3 days before the submission deadline. Responses to request for clarification will be communicated directly in the portal.

Documents to be submitted	<p> <input checked="" type="checkbox"/> Annex 2 Quotation Submission Form duly completed and signed <input checked="" type="checkbox"/> Annex 3 Technical and Financial Offer duly completed and signed and in accordance with the Schedule of Requirements Annex 1, detailing product description / brochures, price quotation per item. <input checked="" type="checkbox"/> Annex 4 Technical Responsiveness Table duly completed and signed <input checked="" type="checkbox"/> Company Profile (short info up to 3 pages), including portfolio/previous corporate experience in similar fields related to the assignment <input checked="" type="checkbox"/> Copy of Company's Registration certificate/articles of incorporation; <input checked="" type="checkbox"/> List and value of minimum 2 (two) similar contracts for the provision and supply of forensic investigation / crime scene kits or similar goods, ongoing or implemented over the last 3 (three) years (copies of the contracts/purchase orders with specification of goods/equipment to be provided): <ul style="list-style-type: none"> • Name of previous contracts • Client & Reference contact details, including active e-mail • Contract Value • Period of activity / Delivery period • Types of goods / equipment delivered </p> <p> <input checked="" type="checkbox"/> Official Letter of Appointment as local representative, if Bidder is submitting a Bid on behalf of an entity located outside the country <input checked="" type="checkbox"/> Manufacturer's Authorization of the Company as a Sales Agent (if Supplier is not the manufacturer) (if applicable); <input checked="" type="checkbox"/> Valid quality Certificate(s) (e.g. ISO 9001, etc., or other quality certification related to scope of the assignment) (if applicable); <input checked="" type="checkbox"/> Detailed description of the offered goods, including design and Technical Data Sheet (including photos, leaflets, brochures and other available materials); <input checked="" type="checkbox"/> Delivery schedule, including delivery time and production time ; <input checked="" type="checkbox"/> Statements or Certificates / Declaration of origin for the offered goods, materials and accessories (where applicable); <input checked="" type="checkbox"/> Statement of warranty as per schedule of requirements; <input checked="" type="checkbox"/> Evidence/Certification of Environmental Sustainability ("Green" Standards) of the Company or the Product being supplied (if any); <input checked="" type="checkbox"/> Certificates / Declarations of Conformity, if any; <input checked="" type="checkbox"/> Statement of satisfactory performance from the top 2 (two) clients in terms of Contract value on delivery of similar goods <input checked="" type="checkbox"/> Financial Statements (Income Statements and Balance Sheets) for the past 2 (two) years (2022 and 2021 or 2023 and 2022, if available) </p>
Evaluation method	<p> <input checked="" type="checkbox"/> The Contract or Purchase Order will be awarded to the lowest price substantially compliant offer </p>
Evaluation criteria	<p> <input checked="" type="checkbox"/> Full compliance with all requirements as specified in Annex 1 <input checked="" type="checkbox"/> Be a legally registered entity or a consortium of firms/organizations <input checked="" type="checkbox"/> Minimum 3 (three) year experience in supplying forensic investigation equipment and consumables; <input checked="" type="checkbox"/> Minimum 2 (two) similar contracts for the provision and supply of forensic investigation / crime scene kits or similar goods, ongoing or implemented over the last 3 (three) years (copies of the contracts/purchase orders with specification of goods/equipment to be provided): <input checked="" type="checkbox"/> Full acceptance of the PO/Contract General Terms and Conditions; <input checked="" type="checkbox"/> Warranty - minimum 12 months per kit including each item included in its composition, validity of consumables not smaller than 12 months. <input checked="" type="checkbox"/> Maximum delivery period not to exceed 45 calendar days upon signature of contract </p>
Type of Contract to be awarded	<p> <input checked="" type="checkbox"/> Contract Face Sheet (Goods and-or Services) </p>

Expected date for contract award	20 May 2024
---	-------------

ANNEX 1 SCHEDULE OF REQUIREMENTS

Technical Specifications: Crime Scene Kits

Item	GENERIC SPECIFICATIONS		UOM	Qty
1.	Crime Scene Marking Kit		kits	4
	Purpose:			
		Used by forensic investigators to accurately document and mark evidence at a crime scene.		
	Composition:			
		1. Case 2. Fiberglass reel tape 3. Metal measuring tape 4. Folding ruler 5. Laser rangefinder		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			
	Case: 1 pc per kit	- ABS - Preferably black or any dark color - Anti corrosive hinges - Two latches - Preferably folding front handle - Supports at least 10 kg - Padding: foam insert or elastic straps		
	Fiberglass reel tape: 1 pc per kit	- Yellow fiberglass; - Length - min. 30 m; - Black millimeter and centimeter markings; - Flexible, resistant to moisture, non-electrically conductive, resistant to physical actions of twisting, bending and breaking; - Versatile end hook; - Rewind crank.		
Metal measuring tape: 1 pc per kit	- Anticorrosive finish; - Yellow color; - Length – 5 m; - Black millimeter and centimeter markings; - Flexible, resistant to moisture, non-electrically conductive, resistant to physical actions of twisting, bending and breaking; - Versatile end hook - Locking mechanism			

Folding ruler: 1 pc per kit	<ul style="list-style-type: none"> - Wooden or metal - Alternative black/white metric marking - Length – 1 m; - Joints with stops 		
Laser rangefinder: 1 pc per kit	<ul style="list-style-type: none"> - Laser beam, class II - Optical device to follow laser beam on long distances - Measuring range: 0.05 m up to 200 m - Accuracy: +/- 2 mm - Backlit display - Tripod mounting possibility - Metric units - Possibility to memorize the last 10 measurements - Continuous measurements function 		
A-Frame evidence marking tents: 1 pc per kit	<ul style="list-style-type: none"> - Letters: A to G - Numbers: 1 to 20 - Yellow color with black printing - Square or rectangular shape of sides - Dimensions: 7-10 cm - Carrying case/bag 		
Evidence marking flags: 1 pc per kit	<ul style="list-style-type: none"> - Metallic pole - Height: 30-40 cm - Colour: red (10 pcs), yellow (10 pcs), green (10 pcs) - Carrying case/bag 		
Photo evidence tape: 3 pcs per kit	<ul style="list-style-type: none"> - Adhesive tape - Width: 10-15 mm - Length: 3 m with 0-30 cm repeat - Metric (millimeters and centimeters) black marking on white matte finish 		
Vinyl magnetic evidence scale: 4 pcs of each color	<ul style="list-style-type: none"> - Metric (centimeters and millimeters) marking - Black on white - White on black - Length – 15 cm 		
Reflective chalk: 3 pcs per kit			
Permanent markers and water-based markers (red, black, blue): 1 pc of each type and color per one kit			
Clutch/mechanical pencil with red, black, blue leads 1 set of leads of each color)			
Fluorescent pencils for forensic marking (red) 3 pcs per kit			
Stickers (set of min 50 pcs): 1 set of each color per kit <ul style="list-style-type: none"> - Red and green - Width – 10 cm - Length – 1.5 cm 			
Documentation:			
	User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		

2.	Universal Forensic Kit		kits	10
	Purpose:			
		A comprehensive set of tools and materials designed for forensic analysis across various types of evidence.		
	Composition:			
		<div>1. Case</div> <div>2. Fingerprint ink pad</div> <div>3. Fingerprint ink roller</div> <div>4. Black fingerprint ink</div> <div>5. Ink cleaning solution or wipes</div> <div>6. Gelatin lifters for fingerprints</div> <div>7. Fingerprint lifters (size 10 x 15 cm)</div> <div>8. Fingerprint lifters (size 5 x 10 cm)</div> <div>9. Transparent lifting tape</div> <div>10. Small particle reagents</div> <div>11. Latent print powders</div> <div>12. Magnetic latent print powders</div> <div>13. Standard fiberglass brush</div> <div>14. Magnetic powder applicator with protective container</div> <div>15. Squirrel or Camel hair brush (type 1)</div> <div>16. Squirrel or Camel hair brush (type 2)</div> <div>17. Digital caliper</div> <div>18. Stainless steel scissors:</div> <div>19. Silicone casting kit</div> <div>20. Silicone rubber spray</div> <div>21. Modeling clay (monochrome bars)</div> <div>22. Round magnet with handle</div> <div>23. Flat tip tweezers (type 1)</div> <div>24. Flat tip tweezers (type 2)</div> <div>25. Sharp tip tweezers</div> <div>26. Surgical scalpel handle</div> <div>27. Disposable sterile scalpel blades</div> <div>28. Metal measuring tape</div> <div>29. Thin sterile surgical gloves (size XL)</div> <div>30. Thick rubber gloves</div> <div>31. Zip-lock plastic bags</div> <div>32. Sterile plastic bags</div> <div>33. Transparent evidence collection plastic tubes with cap</div> <div>34. Stickers with figures 0 to 9</div> <div>35. Stickers with printed arrows</div> <div>36. Universal plier</div> <div>37. Fixing spray</div> <div>38. Plastic mixing container/bowl</div> <div>39. Permanent red marker</div> <div>40. Flashlight</div>		

		41. Magnifying glass 42. Phillips screwdriver 43. Flat-head screwdriver 44. Voltage tester screwdriver 45. Glass cutter 46. Reflective chalk		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			
	Case:	- ABS - Preferably black or any dark color - Anti corrosive hinges - Two latches - Preferably folding front handle - Supports at least 10 kg; - Padding: foam insert or elastic straps		
	Fingerprint ink pad:	- Plastic or similar material case - Snap lock (or similar) lid - Minimum active surface – 12 cm x 6 cm		
	Fingerprint ink roller:	- Made of rubber - Width: 5-7 cm - Handle bracket		
	Black fingerprint ink (min 100 ml)			
	Ink cleaning solution (min 100 ml) or wipes (min 100 pcs)			
	Gelatin lifters for fingerprints: 30 pcs of each color per kit	Size: 13 x 18 cm (± 2 cm) Color: - Black - White - Transparent - Supplied within a cold bag Note: Only 20 pcs of lifters of each color will be supplied within the main case. The remaining quantity of lifters will be supplied in a separate resupply box within a cold bag.		
	Fingerprint lifters: 20 pcs of each color per kit	Size: 10 x 15 cm (± 2 cm) - Black - White - Transparent		
	Fingerprint lifters: 50 pcs of each color per kit	Size: 5 x 10 cm (± 2 cm) - Black - White - Transparent Note: Only 20 pcs of lifters of each color will be		

		supplied within the main case. The remaining quantity of lifters will be supplied in a separate resupply box within a cold bag.		
	Transparent lifting tape:	- Length: min 10 m - Width: min 3 cm, 5 cm, 10 cm		
	Small particle reagents:	- Light-coloured surfaces – 500 ml - Dark-coloured surfaces – 500 ml Each container shall be provided with a dispenser spray head		
	Latent print powders:	- Color: white, black and silver - Application by powdering or by using special brushes - Operating temperature range: -10° C/ +45°C - Storage box		
	Magnetic latent print powders:	- Black - Double contrast - Fluorescent red - Storage box		
	Standard fiberglass brush: 2 pcs per kit	- Length of active part of brush – at least 5 cm - Protective container and cap		
	Magnetic powder applicator with protective container			
	Squirrel or Camel hair brush: 2 pcs per kit	- Uncut hair - Semicircle - Wooden handle – 10-15 cm long - Length of active part of brush – at least 3 cm - Protective container and cap		
	Squirrel or Camel hair brush (type 2): 2 pcs per kit	- Uncut hair - Semicircle - Wooden handle – 15-20 cm long - Length of active part of brush – at least 3 cm - Protective container and cap		
	Digital caliper:	- Stainless steel - Length – 15 cm		
	Stainless steel scissors:	- Rounded tips - Length – 10-15 cm		
	Silicone casting kit:1 set	- Liquid silicone rubber – 1 tube of min 150 g - 1 vial of catalyst – quantity sufficient for the liquid silicone rubber quantity - Plastic spatula – 1 pc - Silicone compound to be of light color, waterproof and able to return to original shape after deformation - Silicone compound hardens within 10 minutes		
	Silicone rubber spray (1 tube of min 150 ml)			
	Modeling clay (monochrome bars, at least 30 g/bar) 3 pcs per kit			
	Round magnet with handle (diameter: 20-25 mm)			
	Flat tip tweezers (type 1):	- Material: Stainless steel		

		- Length: 10-15 cm		
Flat tip tweezers (type 2):		- Material: Stainless steel - Length: 20-25 cm		
Sharp tip tweezers:		- Material: Stainless steel - Length: 8-10 cm		
Surgical scalpel handle: plastic or similar material.				
Disposable sterile scalpel blades (packed separately), compatible with the surgical scalpel handle. 10 pcs per 1 kit				
Metal measuring tape:		- Anticorrosive finish - Yellow color - Length: 5 m - Black millimeter and centimeter markings - Versatile end hook - Locking mechanism		
Thin sterile surgical gloves (size XL) 5 pairs per kit				
Thick rubber gloves				
Zip-lock plastic bags: 5 pcs per kit		- Dimensions: 15-20 cm x 45-50 cm		
Sterile plastic bags: 5 pcs per kit		- Dimensions: 10-15 cm x 40-45 cm		
Transparent evidence collection plastic tubes with cap: 5 pcs per kit		- Diameter: 2 cm - Length: 15-20 cm		
Stickers with figures 0 to 9 (6 pcs per A4 sheet, one sheet for each figure) 10 sheets per kit				
Stickers with printed arrows: 4 sheets		- Black marking on white background - Length: 30-40 mm - Width: 10-20 mm		
Universal plier:		- Anticorrosive finish - Insulated handles		
Fixing spray:		- 1 tube of 100 ml - Colourless - Operation temperature range: -10° C/ +40°C - Drying/hardening time: up to 5 min		
Plastic mixing container/bowl, diameter at the top 12-15 cm				
Permanent red marker				
Flashlight		- UV LED technology - Wave length: 395 nm (+/- 5 nm) - Dimensions: diameter – up to 4 cm, length – up to 15 cm - Rechargeable battery - Battery life (fully charged): at least 2 hours - Battery charger included - Body: metal or highly resistant plastic		

	Magnifying glass:	- Diameter: 75 mm - Magnification: 3.5 x or more		
	Phillips screwdriver: length 15-17 cm			
	Flat-head screwdriver: length 18-20 cm			
	Voltage tester screwdriver: Voltage: AC 100-500V			
	Glass cutter			
	Reflective chalk: 3 pcs per kit			
	Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		
3.	Trace evidence kit		kits	4
	Purpose:			
		Used to collect and preserve small, often microscopic, pieces of physical evidence found at a crime scene.		
	Composition:			
		1. Case 2. MIKROSIL Kit (brown and white) 3. Dental gypsum or Crownstone plaster 4. Rubber mixing container/bowl 5. Adjustable casting frame 6. Wax spray 7. Fixing spray 8. Pig hairbrush 10. Plastic pipette with rubber bulb 11. Filter paper 12. Plastic spatula 13. Metal/wood “junior” hacksaw 14. Plumber wrench 15. Hammer 16. Diagonal 17. Carpenter’s chisel 18. Metal cutting scissors 19. Metal tweezers 20. Metal Kocher forceps 21. Silicone spray 22. Universal Tool (Breaker) - 1 Hole		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			

	Case:	<ul style="list-style-type: none"> - ABS - Preferably black or any dark color - Anti corrosive hinges - Two latches - Preferably folding front handle - Supports at least 10 kg; - Padding: foam insert or elastic straps 		
	MIKROSIL Kit (brown and white):	<p>Material:</p> <ul style="list-style-type: none"> - Plastic or light metal tubes with a capacity between 20-50 ml. <p>Closure System:</p> <ul style="list-style-type: none"> - Tight closing system (caps) to securely hold the contents and prevent leakage during transport and storage. <p>Labeling:</p> <ul style="list-style-type: none"> - Product Name: Clearly labeled with the name of the product. <p>Supplier Information:</p> <ul style="list-style-type: none"> - Name and contact details of the supplier. <p>Quantity:</p> <ul style="list-style-type: none"> - Amount of substance contained in the kit. <p>Expiration Date:</p> <ul style="list-style-type: none"> - Warranty period indicating the expiration date for the product. <p>Safety Information:</p> <ul style="list-style-type: none"> - Mentions regarding the nature of the substance and its toxicity level for the human body. <p>Storage, Transport, and Use Instructions:</p> <ul style="list-style-type: none"> - Clear instructions regarding mandatory conditions for storage, transport, and use to ensure safe handling and optimal performance. 		
	Dental gypsum or Crownstone plaster:	<ul style="list-style-type: none"> - Powder - 1kg bag - Operating temperature range: -10° C/ +45°C - Hardens within 30 min 		
	Rubber mixing container/bowl:	<ul style="list-style-type: none"> - Diameter at the top 12-15 cm - Resistant to chemicals and petroleum based oils 		
	Adjustable casting frame: 3 pcs per kit	<ul style="list-style-type: none"> - Made of lightweight metal or highly resistant plastic - Height of frame: 2-6 cm - Rectangular shape with a possibility for extension on the long side 		

		<ul style="list-style-type: none"> - Size 1 (WxL): 15-20 cm x 30-35 cm (extendable up to 45-50 cm) - Size 2 (WxL): 30-35 cm x 45-50 (extendable up to 60-65 cm) 		
	Wax spray: 3 pcs per kit	<ul style="list-style-type: none"> - Volume 500 ml - Colorless - Operation temperature range: -10° C/ +40°C 		
	Fixing spray: 3 pcs per kit	<ul style="list-style-type: none"> - Volume 500 ml - Colorless - Operation temperature range: -10° C/ +40°C - Drying/hardening time: up to 5 min 		
	Pig hairbrush:	<ul style="list-style-type: none"> - Wooden handle - Length:10-15 cm - Brush length: 3-5 cm 		
	Plastic pipette with rubber bulb			
	Filter paper: : 5 pcs per kit			
	Plastic spatula			
	Metal/wood "junior" hacksaw			
	Plumber wrench			
	Hammer (250 gr) with wooden handle, min 20 cm long			
	Diagonal cutting pliers (min 15 cm long)			
	Carpenter's chisel (min 10 cm long, active part 3 cm)			
	Metal cutting scissors (snips) (20 cm long)			
	Metal tweezers	<ul style="list-style-type: none"> - Length: 10-15 cm Type: <ul style="list-style-type: none"> - Straight edge - Curved point 		
	Kocher forceps	<ul style="list-style-type: none"> -Material: stainless steel Type: straight <ul style="list-style-type: none"> - Length: 10-15 cm 		
	Universal Tool (Breaker) - 1 Hole	Material: <ul style="list-style-type: none"> - High-grade stainless steel or durable alloy construction for strength and longevity. - Corrosion-resistant coating to prevent degradation over time. Design: <ul style="list-style-type: none"> - Ergonomic handle design for comfortable grip and precise control. - Single hole for efficient breaking and separation of materials. - Compact and lightweight for easy handling and portability. Functionality: <ul style="list-style-type: none"> - Capable of breaking a wide range of materials commonly encountered in forensic trace evidence analysis, including glass, fibers, and 		

		<div>ceramics.</div> <div><div>- Precision-engineered hole for controlled braking without excessive fragmentation.</div></div> <div>Compatibility:</div> <div><div>- Compatible with standard forensic laboratory procedures and equipment.</div><div>- Designed to integrate seamlessly into trace evidence kits for comprehensive forensic analysis.</div></div> <div>Durability:</div> <div><div>- Engineered for durability and reliability in demanding forensic environments.</div><div>- Resistant to wear and deformation during repeated use.</div></div> <div>Safety Features:</div> <div><div>- Smooth edges and surfaces to minimize the risk of accidental injury during operation.</div><div>- Non-slip grip for enhanced control and safety.</div></div> <div>Maintenance:</div> <div><div>- Low-maintenance design for ease of upkeep.</div><div>- Recommended periodic cleaning and inspection to ensure optimal performance.</div></div>		
	Silicone spray: 3 pcs per kit	<div>- Volume: 500 ml</div> <div>- Type: colorless</div>		
	Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		
4.	Microparticle evidence kit		kits	3
	Purpose:			
		Used to collect and preserve tiny particles, often invisible to the naked eye, found at a crime scene.		
	Composition:			
		<div>1. Case</div> <div>2. Transparent lifting tape</div> <div>3. Waxed paper</div> <div>4. Portable vacuum cleaner</div>		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			
	Case:	<div>- ABS</div> <div>- Preferably black or any dark color</div> <div>- Anti corrosive hinges</div> <div>- Two latches</div> <div>- Preferably folding front handle</div>		

		<ul style="list-style-type: none">- Supports at least 10 kg;- Padding: foam insert or elastic straps		
	Transparent lifting tape: 1 pc of each dimension per kit	<ul style="list-style-type: none">- Length: 10 m- Width: min 3 cm, 5 cm, 10 cm		
	Waxed paper (A4), min 20 sheets per set, 3 sets per kit			
	Portable vacuum cleaner:	<ul style="list-style-type: none">- Rechargeable battery- Battery life (fully charged) – min 2 hours- Power supply: 210-240 V AC- Filter assembly- Disposable polyester micro-filters – 100 pcs- Lockable plastic case		
5.	Biological traces/particles kit		kits	4
	Purpose:			
		Used to collect and preserve biological materials found at a crime scene, such as blood, saliva, semen, or skin cells		
	Composition:			
		<ol style="list-style-type: none">1. Case2. Unpowdered sterile medical gloves3. Forensic sterile cotton swab with plastic tip protector4. Cardboard box5. DNA controlled absorbent cotton tipped swab6. Surgical scalpel handle (metal)7. Disposable blades for the scalpel8. Sterile syringes with needles9. Sterile gauze pads10. Sterile water11. Solution for sterilization of instruments and container12. Evidence enveloped with printed text13. Paper bags14. Writable stickers15. Sterile cotton swabs16. Disposable plastic tweezers		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			
	Case:	<ul style="list-style-type: none">- ABS- Preferably black or any dark color- Anti corrosive hinges- Two latches		

		<ul style="list-style-type: none">- Preferably folding front handle- Supports at least 10 kg;- Padding: foam insert or elastic straps		
	Unpowdered sterile medical gloves (sizes M and L): 10 pairs of each size per kit			
	Forensic sterile cotton swab with plastic tip protector: 10 pcs of per kit			
	Cardboard box, adapted to the dimensions of the swab to individually store and transport the collected samples : 10 pcs of per kit			
	DNA controlled absorbent cotton tipped swab: 50 pcs per kit			
	Surgical scalpel handle (metal)			
	Disposable sterile scalpel blades (packed separately) compatible with surgical scalpel handle: 10 pcs per kit			
	Sterile syringes with needles, 5 ml: 15 pcs per kit			
	Sterile gauze pads (10x10 cm) (5 pcs per pack)			
	Sterile water (10 ml) 1 bottle per kit			
	Solution for sterilization of instruments (100 ml) and container	<ul style="list-style-type: none">- 90 ml of 10% sodium hypochlorite solution- 10 ml of 98% ethanol		
	Evidence envelopes with printed text: 20 pcs of each dimension per kit	<ul style="list-style-type: none">- White transparent paper- Dimensions:<ul style="list-style-type: none">• 9 x 17 cm• 16 x 23 cm• 10 x 15 cm		
	Paper bags (80 x 100 cm) with printed text: 10 pcs per kit			
	Writable stickers: 20 pcs per kit			
	Sterile cotton swabs: 50 pcs per kit			
	Disposable plastic tweezers: 10 pcs			
	Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		
6.	Electrostatic Dust Print Lifter Kit		kit	1
	Purpose:			
		Used to lift and preserve dust impressions, footprints, or other latent prints from surfaces at a crime scene using an electrostatic charge.		
	Composition:			
		<ul style="list-style-type: none">- Electrostatic Voltage Control- Ground Plane, nickel-plated steel- Ground Plane Polycarbonate Insulating Sheet- Insulated Roller- 9V Alkaline Battery- Metalized Lifting Mats in Protective Tube- Lifting Mats- Lifting Material		

		<div>- Static Discharge Cable - Carrying case</div>		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			
	Carrying case			
	Electrostatic Voltage Control			
	Ground Plane, nickel-plated steel (10 cm x 15 cm) ± 2 cm			
	Ground Plane Polycarbonate Insulating Sheet (12 cm x 17 cm) ± 2 cm			
	Insulated roller			
	9V alkaline battery 2 pcs			
	Metalized lifting mats in protective tube: set of 5 pcs			
	Lifting mats: 2 sets of 10 pcs each			
	Lifting material roll			
	Static discharge cable			
	Carrying case			
	Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		
7.	Crime Scene Sketch Kit		kits	4
	Purpose:			
		Used to accurately document the layout, dimensions, and spatial relationships of a crime scene.		
	Composition:			
		<div>1. Pad of Graph Paper (50 Sheets) 2. Plastic ruler (30 cm) 3. Clutch/mechanical sketch pencil 4. Traffic Template 5. Lavatory Planning Template 6. House Furnishings Template 7. Store Layout Template 8. Office Plan Template 9. Universal Compass 10. Metal measuring tape (15 m) 11. Set of Professional Drawing Instruments 12. Triangle ruler (30°/60°/90°) (transparent plastic) 13. Triangle ruler (45°/45°/90°) (transparent plastic) 14. Sketching and Detail Drawing of the Crime Scene Booklet 15. Carrying case with handle</div>		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered		

		specification.		
	Specifications:			
	Pad of Graph Paper (50 Sheets)			
	Plastic ruler (30 cm)			
	Clutch/mechanical sketch pencil			
	Traffic Template			
	Lavatory Planning Template			
	House Furnishings Template			
	Store Layout Template			
	Office Plan Template			
	Universal Compass			
	Metal measuring tape (15 m)			
	Set of Professional Drawing Instruments			
	Triangle ruler (30°/60°/90°) (transparent plastic)			
	Triangle ruler (45°/45°/90°) (transparent plastic)			
	Sketching and Detail Drawing of the Crime Scene Booklet			
	Carrying case with handle			
	Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		
8.	Portable Fingerprinting Kit		kits	2
	Purpose:			
		Designed to enable law enforcement personnel or forensic investigators to collect fingerprint evidence quickly and efficiently at crime scenes or during investigations.		
	Composition:			
		1. Master carrying case 2. Folding fingerprinting station with a single cardholder and glass slab affixed 3. Tube of grade A fingerprint slab ink 4. Roller for applying ink evenly 5. Cleaning solution for slab and roller 6. Fingerprint ridge builder pump spray 7. Ink remover towelettes 8. Criminal fingerprint record cards 9. Re-print tabs 10. Nitrile barrier gloves		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			
	Master Carrying Case:	A sturdy case for transporting and storing the kit.		

	Folding Fingerprinting Station:	A portable workstation with a cardholder and glass slab for fingerprinting.		
	Grade A Fingerprint Slab Ink:	High-quality ink for fingerprinting.		
	Roller for applying ink evenly:	A roller used for applying fingerprint ink evenly.		
	Cleaning solution for slab and roller:	Cleaner for removing ink residue from the slab and roller.		
	Spray for enhancing 1 bottle per kit fingerprint visibility:	Spray to enhance fingerprint visibility and ridge detail.		
	Ink remover towelettes: 1package per kit	Towelettes for cleaning ink from fingers or surfaces.		
	Criminal fingerprint records:	Cards for recording and preserving fingerprint evidence.		
	Re-Print Tabs:	Tabs used for re-inking fingerprints if necessary.		
	Nitrile Barrier Gloves: 20 pcs per kit	Gloves to prevent contamination and protect personnel during fingerprinting procedures.		
	Documentation			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		
9.	Relief Kit. Papillation and Fingerprint Removal		kits	9
	Purpose:			
		Designed for papillation and fingerprint removal, serves the purpose of aiding forensic investigators in the delicate process of collecting and preserving fingerprint evidence from surfaces.		
	Composition:			
		1. Fingerprint powder 2. Magnetic Latent Fingerprint Powder		
	Accessories:			
		Any other accessory which is required for operation should be included and described in your offered specification.		
	Specifications:			
	Fingerprint powder	- Adhesion to papillary ridge secretions. - Non-adhesion to areas adjacent to papillary ridges. - Homogeneous, non-hygroscopic, and non-agglomerating particles. - Black and fluorescent yellow for visibility under UV rays. - Easy removal from surfaces using simple procedures. - Consistency of physico-chemical characteristics in the temperature range of -10°C to +45°C.		
	Magnetic Latent	Types:		

Fingerprint Powder 1 container for each type per kit	- black double contrast - fluorescent orange - fluorescent green and fluorescent red Weight: - 200 gr per container		
Ostrich/Marabu Brush	Type: - Transparent plastic anti-static sprayer, individually packed. Working length of the brush: - length of at least 80 mm and a thickness at the - thickness at the handle joint of at least 8 mm		
Squirrel/Camel Hair Brushes Minimum Quantity: 3 pieces per set	Shape: Semi-Round Trim: Untrimmed Length: 14 cm and 19 cm (respectively for the two sizes) Handle Material: Wooden Handle Length: Standard Active Part Length: 30-40 mm Active Part Thickness: At least 8 mm at the joint with the handle Packaging: Individual packaging for each brush set		
Glass fiber wire brush: Minimum 4 pieces per set	Active Part Length: Minimum 50 mm Active Part Thickness: Minimum 8 mm at the joint with the handle Handle Material: Plastic or Wooden Handle Length: Standard		
Black anodised aluminum magnetic brush 2 pcs per kit	Equipped with spring rod for demagnetization.		
Master kit small particle reagent (SPR): 2 containers with pump, 250-500 ml	Description: Set of reagents for revealing latent papillary marks on wet, white and black surfaces.		
Gelatine Lifters/ Gelatine sheets for the lift of latent papillary mark: min. 10 sheets of each type per kit	Dimensions: 13x18 cm (± 2 cm) Type: white, black and transparent Components: <u>Support:</u> <ul style="list-style-type: none"> ● Material: Plastic or paper laminated on both sides ● Properties: Resistant to mechanical stretching, malleable, unaffected by humidity conditions, opaque (for black foil) and transparent (for transparent foil), allows writing without affecting adhesive <u>Adhesive Component:</u> <ul style="list-style-type: none"> ● Homogeneous adhesive substance, uniformly distributed over the entire surface 		

		<ul style="list-style-type: none"> ● Adhesive to powders used to reveal papillary marks and protective film ● Adhesion to support without absorption or diffusion of powders ● Stable physical-chemical characteristics (-10°C to +45°C), black (for black foil) and transparent (for transparent foil) ● Allows gluing, removal, and regluing of protective film <u>without affecting adhesive</u> <p><u>Protective Layer:</u></p> <ul style="list-style-type: none"> ● Material: Transparent plastic foil, malleable, no adhesion to embossing powder ● Easily removable from the adhesive surface. 		
	Instant adhesive sheets for lifting latent papillary mark: min. 10 sheets for each type per ki	Dimensions: 10x15 cm (±2 cm) Type: White, black and transparent.		
	Latent papillar mark lift foil set. transfer film and fixation of traces: min 50 sheets for each type per kit	Dimensions: 5x10 cm (±2 cm) Type: White, black and transparent.		
	Adhesive film with mm and cm measuring unit: 2 pcs, each 20 m of each type per kit.	Type: White, black and transparent. Width: 10 - 14 mm		
	Case:	Metal cover with tight magnetic closure		
	Rubber roller:	Width of the active surface: 5-7 cm Handle: Wood		
	Fingerprint Ink 1container per kit Pad:	Solution/gel for cleaning mascara, 1 container of at least 100 ml;		
	Hemident reagent 1 bottle per kit			
	Stainless steel scissors			
	Plastic container for brushes and powder:	Durable plastic material to withstand various environmental conditions and dust		
	Flat tweezers:	Material: stainless steel Length: 23-26 cm		
	Sharp tweezers:	Material: stainless steel Length: 8-10 cm		
	Device for cutting glass:			
	Magnifying glass :	Diameter: 75 mm Magnification: not less than 3.5x		
	Metal measuring tape:	Length: 3m		
	Ink cleaning solution:1 bottle per kit	Volume: 250 ml or bigger		
	Flashlight with clamp:	Illumination:		

		<ul style="list-style-type: none"> - High-intensity light source, providing bright illumination for various applications - Equipped with efficient LED bulbs <p>Clamp Mechanism:</p> <ul style="list-style-type: none"> - The flashlight incorporates a sturdy clamp mechanism for secure attachment to surfaces or objects. - The clamp is designed to provide firm grip and stability, allowing for hands-free operation in diverse environments. - Adjustable clamp size to accommodate a wide range of objects, such as poles, handles, or surfaces of varying thickness. <p>Power Source:</p> <ul style="list-style-type: none"> - Powered by rechargeable batteries. - Offers convenient charging options, such as USB charging ports or docking stations, for hassle-free recharging. 		
	Transparent lifting tape:	Minimum length: 15 m Width: 5-9 cm		
	Transparent lifting tape:	Minimum length: 15 m Width: 3-5 cm		
	Ninhydrin Spray: 1 bottle	Volume: 100 ml or more		
	Rhodamine Crystals 6G	Container Type: Glass container with threaded cap Capacity: 25 grams		
	Gentian Violet Crystals:	Container Type: Glass container with threaded lid Capacity: 25 grams		
	1,8-Diazafluoren-9- one bottle per kit	Type: solution Container Type: Glass or plastic container with screw cap Capacity: 120 milliliters		
	Portable lamp with UV rays:	Power Output: Minimum equivalent power of 6 watts. Wavelength: Emitting UV light with a wavelength of 365 nanometers (nm), suitable for specific fluorescence excitation. Dimensions: Approximately 7-10 cm in width, 30-35 cm in height, and 5-8 cm in depth, providing a compact and portable form factor for easy handling		

		<p>and transportation.</p> <p>Handle: Equipped with a handle for comfortable grip and convenient carrying during use.</p> <p>Power Supply: Operates on either 220 volts AC power supply or battery power.</p> <p>Battery: Includes a rechargeable battery for portable operation.</p> <p>Charger: equipped.</p>		
	Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).		

Delivery Requirements

Delivery Requirements	
Delivery date and time	Bidder shall deliver the goods no later than 45 calendar days after Contract signature.
Delivery Terms (INCOTERMS 2020)	DDP Chisinau, Republic of Moldova
Customs clearance (must be linked to INCOTERM)	<input checked="" type="checkbox"/> Supplier/bidder UNDP shall provide a Tax Exemption letter for Customs clearance
Exact Address(es) of Delivery Location(s)	Chisinau municipality, Republic of Moldova, “Stefan cel Mare Academy”, Address No. 21 Gh. Asachi Street
Packing Requirements	Factory packaging
Training on Operations and Maintenance	Not required The supplier shall provide: User manuals for each kit (if applicable), preferable in Romanian language (optional English or Russian).
Warranty Period	Minimum 12 months per kit including each item included in its composition, validity of consumables not smaller than 12 months.
After-sales service and local service support requirements	Replacement of goods with defects or which do not correspond to the required quality with new goods of similar or better quality within acceptable market delivery standards, however not later than 2 (two) months since official request. Existence of Service centre in Moldova is required (not compulsory to be authorised)

ANNEX 2 QUOTATION SUBMISSION FORM

Bidders are requested to complete this form, including the Company Profile and Bidder's Declaration, sign it and return it as part of their quotation along with Annex 3 Technical and Financial Offer. The Bidder shall fill in this form in accordance with the instructions indicated. No alterations to its format shall be permitted and no substitutions shall be accepted.

Name of Bidder	Click or tap here to enter text.	
RFQ reference	RFQ24/02854: Purchase of Crime Scene Investigation Kits for the Police Academy's Forensic Laboratory	Date Click or tap to enter a date.

Company Profile

Item Description	Detail
Legal name of bidder or Lead entity for JVs	Click or tap here to enter text.
Legal Address, City, Country	Click or tap here to enter text.
Website	Click or tap here to enter text.
Year of Registration	Click or tap here to enter text.
Legal structure	Choose an item.
Are you a UNGM registered vendor?	<input type="checkbox"/> Yes <input type="checkbox"/> No If yes, insert UNGM Vendor Number
Quality Assurance Certification (e.g. ISO 9000 or Equivalent) (If yes, provide a Copy of the valid Certificate)	<input type="checkbox"/> Yes <input type="checkbox"/> No
Does your Company hold any accreditation such as ISO 14001 or ISO 14064 or equivalent related to the environment? (If yes, provide a Copy of the valid Certificate)	<input type="checkbox"/> Yes <input type="checkbox"/> No
Does your Company have a written Statement of its Environmental Policy? (If yes, provide a Copy)	<input type="checkbox"/> Yes <input type="checkbox"/> No
Does your organization demonstrate significant commitment to sustainability through some other means, for	<input type="checkbox"/> Yes <input type="checkbox"/> No

example internal company policy documents on women empowerment, renewable energies or membership of trade institutions promoting such issues <i>(If yes, provide a Copy)</i>	
Is your company a member of the UN Global Compact	<input type="checkbox"/> Yes <input type="checkbox"/> No
Bank Information	Bank Name Click or tap here to enter text. Bank Address Click or tap here to enter text. IBAN Click or tap here to enter text. SWIFT/BIC Click or tap here to enter text. Account Currency Click or tap here to enter text. Bank Account Number Click or tap here to enter text.
Previous relevant experience Minimum 2 (two) similar contracts for the provision and supply of forensic investigation / crime scene kits or similar goods, ongoing or implemented over the last 3 (three) years: <i>(copies of provided contracts should be attached to the quotation)</i>	
Name of previous contracts	Client & Reference Contact Details including e-mail Contract Value Period of activity Types of activities undertaken and goods delivered

Bidder's Declaration

Yes	No	
<input type="checkbox"/>	<input type="checkbox"/>	Requirements and Terms and Conditions I/We have read and fully understand the RFQ, including the RFQ Information and Data, Schedule of Requirements, the General Conditions of Contract, and any Special Conditions of Contract. I/we confirm that the Bidder agrees to be bound by them.
<input type="checkbox"/>	<input type="checkbox"/>	I/We confirm that the Bidder has the necessary capacity, capability, and necessary licenses to fully meet or exceed the Requirements and will be available to deliver throughout the relevant Contract period.
<input type="checkbox"/>	<input type="checkbox"/>	Ethics In submitting this Quote I/we warrant that the bidder has not entered into any improper, illegal, collusive or anti-competitive arrangements with any Competitor; has not directly or indirectly approached any representative of the Buyer (other than the Point of Contact) to lobby or solicit information in relation to the RFQ ;has not attempted to influence, or provide any form of personal inducement, reward or benefit to any representative of the Buyer.
<input type="checkbox"/>	<input type="checkbox"/>	I/We confirm to undertake not to engage in proscribed practices, , or any other unethical practice, with the UN or any other party, and to conduct business in a manner that averts any financial, operational, reputational or other undue risk to the UN and we have read the United Nations Supplier Code of Conduct https://www.un.org/Depts/ptd/about-us/un-supplier-code-conduct and acknowledge that it provides the minimum standards expected of suppliers to the UN.
<input type="checkbox"/>	<input type="checkbox"/>	Conflict of interest I/We warrant that the bidder has no actual, potential, or perceived Conflict of Interest in submitting this Quote or entering a Contract to deliver the Requirements. Where a Conflict

Yes	No	
		of Interest arises during the RFQ process the bidder will report it immediately to the Procuring Organisation's Point of Contact.
<input type="checkbox"/>	<input type="checkbox"/>	Prohibitions, Sanctions I/We hereby declare that our firm, ultimate beneficial owners, its affiliates or subsidiaries or employees, including any JV/Consortium members or subcontractors or suppliers for any part of the contract is not under procurement prohibition by the United Nations, including but not limited to prohibitions derived from the Compendium of United Nations Security Council Sanctions Lists and have not been suspended, debarred, sanctioned or otherwise identified as ineligible by any UN Organization or the World Bank Group or any other international Organization.
<input type="checkbox"/>	<input type="checkbox"/>	Bankruptcy I/We have not declared bankruptcy, are not involved in bankruptcy or receivership proceedings, and there is no judgment or pending legal action against them that could impair their operations in the foreseeable future.
<input type="checkbox"/>	<input type="checkbox"/>	Offer Validity Period I/We confirm that this Quote, including the price, remains open for acceptance for the Offer Validity.
<input type="checkbox"/>	<input type="checkbox"/>	I/We understand and recognize that you are not bound to accept any Quotation you receive, and we certify that the goods offered in our Quotation are new and unused.
<input type="checkbox"/>	<input type="checkbox"/>	By signing this declaration, the signatory below represents, warrants and agrees that he/she has been authorised by the Organization/s to make this declaration on its/their behalf.

Signature _____

Name Click or tap here to enter text.

Title Click or tap here to enter text.

Date Click or tap to enter a date.

ANNEX 3 TECHNICAL AND FINANCIAL OFFER - GOODS

Bidders are requested to complete this form, sign it and return it as part of their bid along with Annex 2 Quotation Submission Form. The Bidder shall fill in this form in accordance with the instructions indicated. No alterations to its format shall be permitted and no substitutions shall be accepted.

Name of Bidder	Click or tap here to enter text.	
RFQ reference	RFQ24/02854: Purchase of Crime Scene Investigation Kits for the Police Academy's Forensic Laboratory	Date Click or tap to enter a date.

Bidders must include all costs related to the performance of the services in their quotation price (delivery of all goods and equipment, transport costs, unloading at the delivery site, installation etc.)

CRIME SCENE INVESTINATION KITS FOR THE POLICE ACADEMY'S FORENSIC LABORATORY

		Currency of the Quotation USD US Dollars for international vendors or MDL Moldovan Leu for local (VAT exclusive)				
		INCOTERMS DDP Chisinau				
Item No	Description	UOM	Qty	Unit price Please insert currency	Other costs (transportation, insurance)	Total price Please insert currency
1	Forensic evidence kit for marking	Kit	4			
2	Universal forensic kit	Kit	10			
3	Forensic Trace Research Kit	Kit	4			
4	Forensic microparticle kit	Kit	3			
5	Forensic kit for collecting biological traces	Kit	4			
6	Electrostatic dust extraction source	Kit	1			
7	Crime scene sketch kit	Kit	4			
8	Portable kit for fingerprinting	Kit	2			
9	Kit for detection, lifting of papillary traces and fingerprints	Kit	9			
Total Price						
Total Final and All-inclusive Price						

*In case the delivery and other costs are included in the value of equipment and cannot be reflected as a separate line of expenses, please specify 0 value in "Delivery and other costs" line in the table above, while the price of 0.01 USD is to be indicated in the section "Lines" of Quantum tender system. In such a situation, the exact Contract price will be corrected during the contract negotiation with the winner.

Compliance with Requirements

	You Responses		
	Yes, we will comply	No, we cannot comply	If you cannot comply, pls. indicate counter - offer
Maximum delivery period not to exceed 45 calendar days upon signature of contract	<input type="checkbox"/>	<input type="checkbox"/>	Click or tap here to enter text.
Minimum Technical Specifications	<input type="checkbox"/>	<input type="checkbox"/>	Click or tap here to enter text.
Warranty period - minimum 12 months per kit including each item included in its composition, validity of consumables not smaller than 12 months.	<input type="checkbox"/>	<input type="checkbox"/>	Click or tap here to enter text.
Delivery Term (INCOTERMS) - DDP Chisinau	<input type="checkbox"/>	<input type="checkbox"/>	Click or tap here to enter text.
Validity of Quotation - 90 calendar days	<input type="checkbox"/>	<input type="checkbox"/>	Click or tap here to enter text.
All provisions of the UNDP General terms and Conditions	<input type="checkbox"/>	<input type="checkbox"/>	Click or tap here to enter text.
Payment terms	<input type="checkbox"/>	<input type="checkbox"/>	Click or tap here to enter text.

Other Information

Estimated weight/volume/dimension of the Consignment	Click or tap here to enter text.
Country/ies of Origin (if export licence required this must be submitted if awarded the contract)	Click or tap here to enter text.

I, the undersigned, certify that I am duly authorized to sign this quotation and bind the company below in event that the quotation is accepted.	
Exact name and address of company Company Name Click or tap here to enter text. Address Click or tap here to enter text. Phone No. Click or tap here to enter text. Email Address Click or tap here to enter text.	Authorized Signature _____ _____ Date Click or tap here to enter text. Name Click or tap here to enter text. Functional Title of Authorised Signatory Click or tap here to enter text. Email Address Click or tap here to enter text.

ANNEX 4 TECHNICAL RESPONSIVENESS TABLE

Bidders are requested to complete this form, sign it and return it as part of their bid along with Annex 2 Quotation Submission Form and Annex 3 Technical and Financial Offer. The Bidder shall fill in this form in accordance with the instructions indicated. No alterations to its format shall be permitted and no substitutions shall be accepted.

Name of Bidder	Click or tap here to enter text.	
RFQ reference	RFQ24/02854: Purchase of Crime Scene Investigation Kits for the Police Academy's Forensic Laboratory	Date Click or tap to enter a date.

Bidders shall supplement the information provided in the tables below with detailed description of the offered goods, including design and Technical Data Sheet (including photos, leaflets, brochures and other available materials)

Technical requirements (Features at least) Conformity		(Yes/No)	Offered Brand, Model and Exact characteristics <small>(Technical characteristics should not be less/worse than those listed in Annex 1 Schedule of Requirements)</small>	
1	Crime Scene Marking Kit		Crime Scene Marking Kit	
	Purpose:		Purpose:	
	Used by forensic investigators to accurately document and mark evidence at a crime scene.		Used by forensic investigators to accurately document and mark evidence at a crime scene.	
	Composition:		Composition:	
	Accessories:		Accessories:	
	Any other accessory which is required for operation should be included and described in your offered specification.			
	Specifications:		Specifications:	

	Case: 1 pc per kit	<ul style="list-style-type: none"> - ABS - Preferably black or any dark color - Anti corrosive hinges - Two latches - Preferably folding front handle - Supports at least 10 kg - Padding: foam insert or elastic straps 		Case: 1 pc per kit	
	Fiberglass reel tape: 1 pc per kit	<ul style="list-style-type: none"> - Yellow fiberglass; - Length - min. 30 m; - Black millimeter and centimeter markings; - Flexible, resistant to moisture, non-electrically conductive, resistant to physical actions of twisting, bending and breaking; - Versatile end hook; - Rewind crank. 		Fiberglass reel tape: 1 pc per kit	
	Metal measuring tape: 1 pc per kit	<ul style="list-style-type: none"> - Anticorrosive finish; - Yellow color; - Length – 5 m; - Black millimeter and centimeter markings; - Flexible, resistant to moisture, non-electrically conductive, resistant to physical actions of twisting, bending and breaking; - Versatile end hook - Locking mechanism 		Metal measuring tape: 1 pc per kit	
	Folding ruler: 1 pc per kit	<ul style="list-style-type: none"> - Wooden or metal - Alternative black/white metric marking - Length – 1 m; - Joints with stops 		Folding ruler: 1 pc per kit	
	Laser rangefinder: 1 pc per kit	<ul style="list-style-type: none"> - Laser beam, class II - Optical device to follow laser beam 		Laser rangefinder: 1 pc per kit	

		on long distances - Measuring range: 0.05m up to 200m - Accuracy: +/- 2 mm - Backlit display - Tripod mounting possibility - Metric units - Possibility to memorize the last 10 measurements - Continuous measurements function			
	A-Frame evidence marking tents: 1 pc per kit	- Letters: A to G - Numbers: 1 to 20 - Yellow color with black printing - Square or rectangular shape of sides - Dimensions: 7-10 cm - Carrying case/bag		A-Frame evidence marking tents: 1 pc per kit	
	Evidence marking flags: 1 pc per kit	- Metallic pole - Height: 30-40 cm - Colour: red (10 pcs), yellow (10 pcs), green (10 pcs) - Carrying case/bag		Evidence marking flags: 1 pc per kit	
	Photo evidence tape: 3 pcs per kit	- Adhesive tape - Width: 10-15 mm - Length: 3 m with 0-30 cm repeat - Metric (millimeters and centimeters) black marking on white matte finish		Photo evidence tape: 3 pcs per kit	
	Vinyl magnetic evidence scale: 4 pcs of each color	- Metric (centimeters and millimeters) marking - Black on white - White on black - Length – 15 cm		Vinyl magnetic evidence scale: 4 pcs of each color	
	Reflective chalk: 3 pcs per kit			Reflective chalk: 3 pcs per kit	
	Permanent markers and water-based markers (red, black,			Permanent markers and	

	blue): 1 pc of each type and color per one kit			water-based markers (red, black, blue): 1 pc of each type and color per one kit			
	Clutch/mechanical pencil with red, black, blue leads 1 set of leads of each color)			Clutch/mechanical pencil with red, black, blue leads 1 set of leads of each color)			
	Fluorescent pencils for forensic marking (red) 3 pcs per kit			Fluorescent pencils for forensic marking (red) 3 pcs per kit			
	Stickers (set of min 50 pcs): 1 set of each color per kit - Red and green - Width – 10 cm - Length – 1.5 cm			Stickers (set of min 50 pcs): 1 set of each color per kit - Red and green - Width – 10 cm - Length – 1.5 cm			
	Documentation:			Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).					
2	Universal Forensic Kit			Universal Forensic Kit			
	Purpose:			Purpose:			
		A comprehensive set of tools and materials designed for forensic analysis across various types of evidence.					
	Composition:			Composition:		Composition:	
		1. Case 2. Fingerprint ink pad 3. Fingerprint ink roller 4. Black fingerprint ink 5. Ink cleaning solution or wipes 6. Gelatin lifters for fingerprints 7. Fingerprint lifters (size 10 x 15 cm) 8. Fingerprint lifters (size 5 x 10 cm)					

		9. Transparent lifting tape 10. Small particle reagents 11. Latent print powders 12. Magnetic latent print powders 13. Standard fiberglass brush 14. Magnetic powder applicator with protective container 15. Squirrel or Camel hair brush (type 1) 16. Squirrel or Camel hair brush (type 2) 17. Digital caliper 18. Stainless steel scissors: 19. Silicone casting kit 20. Silicone rubber spray 21. Modeling clay (monochrome bars) 22. Round magnet with handle 23. Flat tip tweezers (type 1) 24. Flat tip tweezers (type 2) 25. Sharp tip tweezers 26. Surgical scalpel handle 27. Disposable sterile scalpel blades 28. Metal measuring tape 29. Thin sterile surgical gloves (size XL) 30. Thick rubber gloves 31. Zip-lock plastic bags 32. Sterile plastic bags 33. Transparent evidence collection plastic tubes with cap 34. Stickers with figures 0 to 9 35. Stickers with printed arrows 36. Universal plier 37. Fixing spray 38. Plastic mixing container/bowl			
--	--	--	--	--	--

	39. Permanent red marker 40. Flashlight 41. Magnifying glass 42. Phillips screwdriver 43. Flat-head screwdriver 44. Voltage tester screwdriver 45. Glass cutter 46. Reflective chalk		
Accessories:		Accessories:	Accessories:
	Any other accessory which is required for operation should be included and described in your offered specification.		
Specifications:		Specifications:	Specifications:
Case:	<ul style="list-style-type: none"> - ABS - Preferably black or any dark color - Anti corrosive hinges - Two latches - Preferably folding front handle - Supports at least 10 kg; - Padding: foam insert or elastic straps 	Case:	
Fingerprint ink pad:	<ul style="list-style-type: none"> - Plastic or similar material case - Snap lock (or similar) lid - Minimum active surface – 12 cm x 6 cm 	Fingerprint ink pad:	
Fingerprint ink roller:	<ul style="list-style-type: none"> - Made of rubber - Width: 5-7 cm - Handle bracket 	Fingerprint ink roller:	
Black fingerprint ink (min 100 ml)		Black fingerprint ink (min 100 ml)	
Ink cleaning solution (min 100 ml) or wipes (min 100 pcs)		Ink cleaning solution (min 100 ml) or wipes (min 100 pcs)	
Gelatin lifters for fingerprints: 30 pcs	Size: 13 x 18 cm (±2 cm) Color:	Gelatin lifters for fingerprints: 30 pcs of each color per kit	

	of each color per kit	<ul style="list-style-type: none"> - Black - White - Transparent - Supplied within a cold bag <p>Note: Only 20 pcs of lifters of each color will be supplied within the main case. The remaining quantity of lifters will be supplied in a separate resupply box within a cold bag.</p>			
	Fingerprint lifters: 20 pcs of each color per kit	Size: 10 x 15 cm (±2 cm) <ul style="list-style-type: none"> - Black - White - Transparent 		Fingerprint lifters: 20 pcs of each color per kit	
	Fingerprint lifters: 50 pcs of each color per kit	Size: 5 x 10 cm (±2 cm) <ul style="list-style-type: none"> - Black - White - Transparent <p>Note: Only 20 pcs of lifters of each color will be supplied within the main case. The remaining quantity of lifters will be supplied in a separate resupply box within a cold bag.</p>		Fingerprint lifters: 50 pcs of each color per kit	
	Transparent lifting tape:	<ul style="list-style-type: none"> - Length: min 10 m - Width: min 3 cm, 5 cm, 10 cm 		Transparent lifting tape:	
	Small particle reagents:	<ul style="list-style-type: none"> - Light-coloured surfaces – 500 ml - Dark-coloured surfaces – 500 ml <p>Each container shall be provided with a dispenser spray head</p>		Small particle reagents:	
	Latent print powders:	<ul style="list-style-type: none"> - Color: white, black and silver - Application by powdering or by using special brushes - Operating temperature range: -10° C/ +45°C - Storage box 		Latent print powders:	

	Magnetic latent print powders:	<ul style="list-style-type: none"> - Black - Double contrast - Fluorescent red - Storage box 		Magnetic latent print powders:	
	Standard fiberglass brush: 2 pcs per kit	<ul style="list-style-type: none"> - Length of active part of brush – at least 5 cm - Protective container and cap 		Standard fiberglass brush: 2 pcs per kit	
	Magnetic powder applicator with protective container			Magnetic powder applicator with protective container	
	Squirrel or Camel hair brush: 2 pcs per kit	<ul style="list-style-type: none"> - Uncut hair - Semicircle - Wooden handle – 10-15 cm long - Length of active part of brush – at least 3 cm - Protective container and cap 		Squirrel or Camel hair brush: 2 pcs per kit	
	Squirrel or Camel hair brush (type 2): 2 pcs per kit	<ul style="list-style-type: none"> - Uncut hair - Semicircle - Wooden handle – 15-20 cm long - Length of active part of brush – at least 3 cm - Protective container and cap 		Squirrel or Camel hair brush (type 2): 2 pcs per kit	
	Digital caliper:	<ul style="list-style-type: none"> - Stainless steel - Length – 15 cm 		Digital caliper:	
	Stainless steel scissors:	<ul style="list-style-type: none"> - Rounded tips - Length – 10-15 cm 		Stainless steel scissors:	
	Silicone casting kit:1 set	<ul style="list-style-type: none"> - Liquid silicone rubber – 1 tube of min 150 g - 1 vial of catalyst – quantity sufficient for the liquid silicone rubber quantity - Plastic spatula – 1 pc - Silicone compound to be of light color, waterproof and able to return to original shape after deformation - Silicone compound hardens within 		Silicone casting kit:1 set	

	10 minutes		
Silicone rubber spray (1 tube of min 150 ml)		Silicone rubber spray (1 tube of min 150 ml)	
Modeling clay (monochrome bars, at least 30 g/bar) 3 pcs per kit		Modeling clay (monochrome bars, at least 30 g/bar) 3 pcs per kit	
Round magnet with handle (diameter: 20-25 mm)		Round magnet with handle (diameter: 20-25 mm)	
Flat tip tweezers (type 1):	- Material: Stainless steel - Length: 10-15 cm	Flat tip tweezers (type 1):	
Flat tip tweezers (type 2):	- Material: Stainless steel - Length: 20-25 cm	Flat tip tweezers (type 2):	
Sharp tip tweezers:	- Material: Stainless steel - Length: 8-10 cm	Sharp tip tweezers:	
Surgical scalpel handle: plastic or similar material.		Surgical scalpel handle: plastic or similar material.	
Disposable sterile scalpel blades (packed separately), compatible with the surgical scalpel handle. 10 pcs per 1 kit		Disposable sterile scalpel blades (packed separately), compatible with the surgical scalpel handle. 10 pcs per 1 kit	
Metal measuring tape:	- Anticorrosive finish - Yellow color - Length: 5 m - Black millimeter and centimeter markings - Versatile end hook - Locking mechanism	Metal measuring tape:	
Thin sterile surgical gloves (size XL) 5 pairs per kit		Thin sterile surgical gloves (size XL) 5 pairs per kit	
Thick rubber gloves		Thick rubber gloves	
Zip-lock plastic bags: 5 pcs per kit	- Dimensions: 15-20 cm x 45-50 cm	Zip-lock plastic bags: 5 pcs per kit	

	Sterile plastic bags: 5 pcs per kit	- Dimensions: 10-15 cm x 40-45 cm		Sterile plastic bags: 5 pcs per kit	
	Transparent evidence collection plastic tubes with cap: 5 pcs per kit	- Diameter: 2 cm - Length: 15-20 cm		Transparent evidence collection plastic tubes with cap: 5 pcs per kit	
	Stickers with figures 0 to 9 (6 pcs per A4 sheet, one sheet for each figure) 10 sheets per kit			Stickers with figures 0 to 9 (6 pcs per A4 sheet, one sheet for each figure) 10 sheets per kit	
	Stickers with printed arrows: 4 sheets	- Black marking on white background - Length: 30-40 mm - Width: 10-20 mm		Stickers with printed arrows: 4 sheets	
	Universal plier:	- Anticorrosive finish - Insulated handles		Universal plier:	
	Fixing spray:	- 1 tube of 100 ml - Colourless - Operation temperature range: -10° C/ +40°C - Drying/hardening time: up to 5 min		Fixing spray:	
	Plastic mixing container/bowl, diameter at the top 12-15 cm			Plastic mixing container/bowl, diameter at the top 12-15 cm	
	Permanent red marker			Permanent red marker	
Flashlight	- UV LED technology - Wave length: 395 nm (+/- 5 nm) - Dimensions: diameter – up to 4 cm, length – up to 15 cm - Rechargeable battery - Battery life (fully charged): at least 2 hours - Battery charger included - Body: metal or highly resistant		Flashlight		

		plastic			
	Magnifying glass:	- Diameter: 75 mm - Magnification: 3.5 x or more		Magnifying glass:	
	Phillips screwdriver: length 15-17 cm			Phillips screwdriver: length 15-17 cm	
	Flat-head screwdriver: length 18-20 cm			Flat-head screwdriver: length 18-20 cm	
	Voltage tester screwdriver: Voltage: AC 100-500V			Voltage tester screwdriver: Voltage: AC 100-500V	
	Glass cutter			Glass cutter	
	Reflective chalk: 3 pcs per kit			Reflective chalk: 3 pcs per kit	
	Documentation:			Documentation:	
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).			
3	Trace evidence kit			Trace evidence kit	
	Purpose:			Purpose:	
		Used to collect and preserve small, often microscopic, pieces of physical evidence found at a crime scene.			
	Composition:			Composition:	
		1. Case 2. MIKROSIL Kit (brown and white) 3. Dental gypsum or Crownstone plaster 4. Rubber mixing container/bowl 5. Adjustable casting frame 6. Wax spray 7. Fixing spray 8. Pig hairbrush 10. Plastic pipette with rubber bulb			

		11. Filter paper 12. Plastic spatula 13. Metal/wood “junior” hacksaw 14. Plumber wrench 15. Hammer 16. Diagonal 17. Carpenter’s chisel 18. Metal cutting scissors 19. Metal tweezers 20. Metal Kocher forceps 21. Silicone spray 22. Universal Tool (Breaker) - 1 Hole			
	Accessories:			Accessories:	
		Any other accessory which is required for operation should be included and described in your offered specification.			
	Specifications:			Specifications:	
	Case:	- ABS - Preferably black or any dark color - Anti corrosive hinges - Two latches - Preferably folding front handle - Supports at least 10 kg; - Padding: foam insert or elastic straps		Case:	
	MIKROSIL Kit (brown and white):	Material: - Plastic or light metal tubes with a capacity between 20-50 ml. Closure System:		MIKROSIL Kit (brown and white):	

		<ul style="list-style-type: none"> - Tight closing system (caps) to securely hold the contents and prevent leakage during transport and storage. <p>Labeling:</p> <ul style="list-style-type: none"> - Product Name: Clearly labeled with the name of the product. <p>Supplier Information:</p> <ul style="list-style-type: none"> - Name and contact details of the supplier. <p>Quantity:</p> <ul style="list-style-type: none"> - Amount of substance contained in the kit. <p>Expiration Date:</p> <ul style="list-style-type: none"> - Warranty period indicating the expiration date for the product. <p>Safety Information:</p> <ul style="list-style-type: none"> - Mentions regarding the nature of the substance and its toxicity level for the human body. <p>Storage, Transport, and Use Instructions:</p> <ul style="list-style-type: none"> - Clear instructions regarding mandatory conditions for storage, transport, and use to ensure safe handling and optimal performance. 			
	Dental gypsum or Crownstone plaster:	<ul style="list-style-type: none"> - Powder - 1kg bag - Operating temperature range: -10° C/ +45°C - Hardens within 30 min 		Dental gypsum or Crownstone plaster:	

	Rubber mixing container/bowl:	<ul style="list-style-type: none"> - Diameter at the top 12-15 cm - Resistant to chemicals and petroleum based oils 		Rubber mixing container/bowl:	
	Adjustable casting frame: 3 pcs per kit	<ul style="list-style-type: none"> - Made of lightweight metal or highly resistant plastic - Height of frame: 2-6 cm - Rectangular shape with a possibility for extension on the long side - Size 1 (WxL): 15-20 cm x 30-35 cm (extendable up to 45-50 cm) - Size 2 (WxL): 30-35 cm x 45-50 (extendable up to 60-65 cm) 		Adjustable casting frame: 3 pcs per kit	
	Wax spray: 3 pcs per kit	<ul style="list-style-type: none"> - Volume 500 ml - Colorless - Operation temperature range: -10° C/ +40°C 		Wax spray: 3 pcs per kit	
	Fixing spray: 3 pcs per kit	<ul style="list-style-type: none"> - Volume 500 ml - Colorless - Operation temperature range: -10° C/ +40°C - Drying/hardening time: up to 5 min 		Fixing spray: 3 pcs per kit	
	Pig hairbrush:	<ul style="list-style-type: none"> - Wooden handle - Length:10-15 cm - Brush length: 3-5 cm 		Pig hairbrush:	
	Plastic pipette with rubber bulb			Plastic pipette with rubber bulb	
	Filter paper: : 5 pcs per kit			Filter paper: : 5 pcs per kit	
	Plastic spatula			Plastic spatula	
	Metal/wood "junior" hacksaw			Metal/wood "junior" hacksaw	
	Plumber wrench			Plumber wrench	
	Hammer (250 gr) with wooden handle, min 20 cm long			Hammer (250 gr) with wooden handle, min 20 cm long	

	Diagonal cutting pliers (min 15 cm long)			Diagonal cutting pliers (min 15 cm long)	
	Carpenter's chisel (min 10 cm long, active part 3 cm)			Carpenter's chisel (min 10 cm long, active part 3 cm)	
	Metal cutting scissors (snips) (20 cm long)			Metal cutting scissors (snips) (20 cm long)	
	Metal tweezers	- Length: 10-15 cm Type: - Straight edge - Curved point		Metal tweezers	
	Kocher forceps	-Material: stainless steel Type: straight - Length: 10-15 cm		Kocher forceps	
	Universal Tool (Breaker) - 1 Hole	Material: - High-grade stainless steel or durable alloy construction for strength and longevity. - Corrosion-resistant coating to prevent degradation over time. Design: - Ergonomic handle design for comfortable grip and precise control. - Single hole for efficient breaking and separation of materials. - Compact and lightweight for easy handling and portability. Functionality: - Capable of breaking a wide range of materials commonly encountered in forensic trace evidence analysis, including glass, fibers, and ceramics.		Universal Tool (Breaker) - 1 Hole	

		<ul style="list-style-type: none"> - Precision-engineered hole for controlled braking without excessive fragmentation. <p>Compatibility:</p> <ul style="list-style-type: none"> - Compatible with standard forensic laboratory procedures and equipment. - Designed to integrate seamlessly into trace evidence kits for comprehensive forensic analysis. <p>Durability:</p> <ul style="list-style-type: none"> - Engineered for durability and reliability in demanding forensic environments. - Resistant to wear and deformation during repeated use. <p>Safety Features:</p> <ul style="list-style-type: none"> - Smooth edges and surfaces to minimize the risk of accidental injury during operation. - Non-slip grip for enhanced control and safety. <p>Maintenance:</p> <ul style="list-style-type: none"> - Low-maintenance design for ease of upkeep. - Recommended periodic cleaning and inspection to ensure optimal performance. 			
	Silicone spray: 3 pcs per kit	- Volume: 500 ml - Type: colorless		Silicone spray: 3 pcs per kit	
		User manuals for the kit (if applicable), preferable in Romanian			

		language (optional English or Russian).				
4	Microparticle evidence kit			Microparticle evidence kit		
	Purpose:			Purpose:		
		Used to collect and preserve tiny particles, often invisible to the naked eye, found at a crime scene.				
	Composition:			Composition:		
		1. Case 2. Transparent lifting tape 3. Waxed paper 4. Portable vacuum cleaner				
	Accessories:			Accessories:		
		Any other accessory which is required for operation should be included and described in your offered specification.				
	Specifications:			Specifications:		
	Case:	- ABS - Preferably black or any dark color - Anti corrosive hinges - Two latches - Preferably folding front handle - Supports at least 10 kg; - Padding: foam insert or elastic straps		Case:		
	Transparent lifting tape: 1 pc of each dimension per kit	- Length: 10 m - Width: min 3 cm, 5 cm, 10 cm		Transparent lifting tape: 1 pc of each dimension per kit		
Waxed paper (A4), min 20 sheets per set, 3 sets per kit		Waxed paper (A4), min 20 sheets per set, 3 sets per kit				
Portable vacuum cleaner:	- Rechargeable battery - Battery life (fully charged) – min 2 hours - Power supply: 210-240 V AC	Portable vacuum cleaner:				

		<ul style="list-style-type: none">- Filter assembly- Disposable polyester micro-filters – 100 pcs- Lockable plastic case			
5	Biological traces/particles kit			Biological traces/particles kit	
	Purpose:			Purpose:	
		Used to collect and preserve biological materials found at a crime scene, such as blood, saliva, semen, or skin cells			
	Composition:			Composition:	
		1. Case 2. Unpowdered sterile medical gloves 3. Forensic sterile cotton swab with plastic tip protector 4. Cardboard box 5. DNA controlled absorbent cotton tipped swab 6. Surgical scalpel handle (metal) 7. Disposable blades for the scalpel 8. Sterile syringes with needles 9. Sterile gauze pads 10. Sterile water 11. Solution for sterilization of instruments and container 12. Evidence enveloped with printed text 13. Paper bags 14. Writable stickers 15. Sterile cotton swabs 16. Disposable plastic tweezers			
	Accessories:			Accessories:	
	Any other accessory which is required for operation should be included and described in your offered specification.				

	Specifications:			Specifications:	
	Case:	<ul style="list-style-type: none">- ABS- Preferably black or any dark color- Anti corrosive hinges- Two latches- Preferably folding front handle- Supports at least 10 kg;- Padding: foam insert or elastic straps		Case:	
	Unpowdered sterile medical gloves (sizes M and L): 10 pairs of each size per kit			Unpowdered sterile medical gloves (sizes M and L): 10 pairs of each size per kit	
	Forensic sterile cotton swab with plastic tip protector: 10 pcs of per kit			Forensic sterile cotton swab with plastic tip protector: 10 pcs of per kit	
	Cardboard box, adapted to the dimensions of the swab to individually store and transport the collected samples : 10 pcs of per kit			Cardboard box, adapted to the dimensions of the swab to individually store and transport the collected samples : 10 pcs of per kit	
	DNA controlled absorbent cotton tipped swab: 50 pcs per kit			DNA controlled absorbent cotton tipped swab: 50 pcs per kit	
	Surgical scalpel handle (metal)			Surgical scalpel handle (metal)	
	Disposable sterile scalpel blades (packed separately) compatible with surgical scalpel handle: 10 pcs per kit			Disposable sterile scalpel blades (packed separately) compatible with surgical scalpel handle: 10 pcs per kit	
	Sterile syringes with needles, 5 ml: 15 pcs per kit			Sterile syringes with needles, 5 ml: 15 pcs per kit	
	Sterile gauze pads (10x10 cm) (5 pcs per pack)			Sterile gauze pads (10x10 cm) (5 pcs per pack)	
	Sterile water (10 ml) 1 bottle per kit			Sterile water (10 ml) 1 bottle per kit	

	Solution for sterilization of instruments (100 ml) and container	- 90 ml of 10% sodium hypochlorite solution - 10 ml of 98% ethanol		Solution for sterilization of instruments (100 ml) and container	
	Evidence envelopes with printed text: 20 pcs of each dimension per kit	- White transparent paper - Dimensions: <ul style="list-style-type: none">• 9 x 17 cm• 16 x 23 cm• 10 x 15 cm		Evidence envelopes with printed text: 20 pcs of each dimension per kit	
	Paper bags (80 x 100 cm) with printed text: 10 pcs per kit			Paper bags (80 x 100 cm) with printed text: 10 pcs per kit	
	Writable stickers: 20 pcs per kit			Writable stickers: 20 pcs per kit	
	Sterile cotton swabs: 50 pcs per kit			Sterile cotton swabs: 50 pcs per kit	
	Disposable plastic tweezers: 10 pcs			Disposable plastic tweezers: 10 pcs	
	Documentation:			Documentation:	
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).			
6	Electrostatic Dust Print Lifter Kit			Electrostatic Dust Print Lifter Kit	
	Purpose:			Purpose:	
		Used to lift and preserve dust impressions, footprints, or other latent prints from surfaces at a crime scene using an electrostatic charge.			
	Composition:			Composition:	
		- Electrostatic Voltage Control - Ground Plane, nickel-plated steel - Ground Plane Polycarbonate Insulating Sheet			

		<ul style="list-style-type: none">- Insulated Roller- 9V Alkaline Battery- Metalized Lifting Mats in Protective Tube- Lifting Mats- Lifting Material- Static Discharge Cable- Carrying case			
	Accessories:			Accessories:	
		Any other accessory which is required for operation should be included and described in your offered specification.			
	Specifications:			Specifications:	
	Carrying case			Carrying case	
	Electrostatic Voltage Control			Electrostatic Voltage Control	
	Ground Plane, nickel-plated steel (10 cm x 15 cm) ± 2 cm			Ground Plane, nickel-plated steel (10 cm x 15 cm) ± 2 cm	
	Ground Plane Polycarbonate Insulating Sheet (12 cm x 17 cm) ± 2 cm			Ground Plane Polycarbonate Insulating Sheet (12 cm x 17 cm) ± 2 cm	
	Insulated roller			Insulated roller	
	9V alkaline battery 2 pcs			9V alkaline battery 2 pcs	
	Metalized lifting mats in protective tube: set of 5 pcs			Metalized lifting mats in protective tube: set of 5 pcs	
	Lifting mats: 2 sets of 10 pcs each			Lifting mats: 2 sets of 10 pcs each	
	Lifting material roll			Lifting material roll	
	Static discharge cable			Static discharge cable	
	Carrying case			Carrying case	
	Documentation:			Documentation:	
		User manuals for the kit (if			

		applicable), preferable in Romanian language (optional English or Russian).			
7	Crime Scene Sketch Kit		kits Purpose: Composition: Accessories: Documentation:	Crime Scene Sketch Kit	
	Purpose:			Purpose:	
		Used to accurately document the layout, dimensions, and spatial relationships of a crime scene.			
	Composition:			Composition:	
		1. Pad of Graph Paper (50 Sheets) 2. Plastic ruler (30 cm) 3. Clutch/mechanical sketch pencil 4. Traffic Template 5. Lavatory Planning Template 6. House Furnishings Template 7. Store Layout Template 8. Office Plan Template 9. Universal Compass 10. Metal measuring tape (15 m) 11. Set of Professional Drawing Instruments 12. Triangle ruler (30°/60°/90°) (transparent plastic) 13. Triangle ruler (45°/45°/90°) (transparent plastic) 14. Sketching and Detail Drawing of the Crime Scene Booklet 15. Carrying case with handle			
	Accessories:			Accessories:	
		Any other accessory which is required for operation should be included and described in your offered specification.			
	Specifications:			Specifications:	

	Pad of Graph Paper (50 Sheets)			Pad of Graph Paper (50 Sheets)		
	Plastic ruler (30 cm)			Plastic ruler (30 cm)		
	Clutch/mechanical sketch pencil			Clutch/mechanical sketch pencil		
	Traffic Template			Traffic Template		
	Lavatory Planning Template			Lavatory Planning Template		
	House Furnishings Template			House Furnishings Template		
	Store Layout Template			Store Layout Template		
	Office Plan Template			Office Plan Template		
	Universal Compass			Universal Compass		
	Metal measuring tape (15 m)			Metal measuring tape (15 m)		
	Set of Professional Drawing Instruments			Set of Professional Drawing Instruments		
	Triangle ruler (30°/60°/90°) (transparent plastic)			Triangle ruler (30°/60°/90°) (transparent plastic)		
	Triangle ruler (45°/45°/90°) (transparent plastic)			Triangle ruler (45°/45°/90°) (transparent plastic)		
	Sketching and Detail Drawing of the Crime Scene Booklet			Sketching and Detail Drawing of the Crime Scene Booklet		
	Carrying case with handle			Carrying case with handle		
	Documentation:		Documentation:			
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).				
8	Portable Fingerprinting Kit			Portable Fingerprinting Kit		
	Purpose:			Purpose:		
		Designed to enable law enforcement personnel or forensic investigators to collect fingerprint evidence quickly				

		and efficiently at crime scenes or during investigations.			
	Composition:			Composition:	
		1. Master carrying case 2. Folding fingerprinting station with a single cardholder and glass slab affixed 3. Tube of grade A fingerprint slab ink 4. Roller for applying ink evenly 5. Cleaning solution for slab and roller 6. Fingerprint ridge builder pump spray 7. Ink remover towelettes 8. Criminal fingerprint record cards 9. Re-print tabs 10. Nitrile barrier gloves			
	Accessories:			Accessories:	
		Any other accessory which is required for operation should be included and described in your offered specification.			
	Specifications:			Specifications:	
	Master Carrying Case:	A sturdy case for transporting and storing the kit.		Master Carrying Case:	
	Folding Fingerprinting Station:	A portable workstation with a cardholder and glass slab for fingerprinting.		Folding Fingerprinting Station:	
	Grade A Fingerprint Slab Ink:	High-quality ink for fingerprinting.		Grade A Fingerprint Slab Ink:	
	Roller for applying ink evenly:	A roller used for applying fingerprint ink evenly.		Roller for applying ink evenly:	
	Cleaning solution for slab and roller:	Cleaner for removing ink residue from the slab and roller.		Cleaning solution for slab and roller:	

	Spray for enhancing 1 bottle per kit fingerprint visibility:	Spray to enhance fingerprint visibility and ridge detail.		Spray for enhancing 1 bottle per kit fingerprint visibility:	
	Ink remover towelettes: 1package per kit	Towelettes for cleaning ink from fingers or surfaces.		Ink remover towelettes: 1package per kit	
	Criminal fingerprint records:	Cards for recording and preserving fingerprint evidence.		Criminal fingerprint records:	
	Re-Print Tabs:	Tabs used for re-inking fingerprints if necessary.		Re-Print Tabs:	
	Nitrile Barrier Gloves: 20 pcs per kit	Gloves to prevent contamination and protect personnel during fingerprinting procedures.		Nitrile Barrier Gloves: 20 pcs per kit	
	Documentation			Documentation	
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).			
9	Relief Kit. Papillation and Fingerprint Removal			Relief Kit. Papillation and Fingerprint Removal	
	Purpose:			Purpose:	
		Designed for papillation and fingerprint removal, serves the purpose of aiding forensic investigators in the delicate process of collecting and preserving fingerprint evidence from surfaces.			
	Composition:			Composition:	
		1. Fingerprint powder 2. Magnetic Latent Fingerprint Powder			
	Accessories:			Accessories:	
	Any other accessory which is required for operation should be included and				

		described in your offered specification.			
	Specifications:			Specifications:	
	Fingerprint powder	- Adhesion to papillary ridge secretions. - Non-adhesion to areas adjacent to papillary ridges. - Homogeneous, non-hygroscopic, and non-agglomerating particles. - Black and fluorescent yellow for visibility under UV rays. - Easy removal from surfaces using simple procedures. - Consistency of physico-chemical characteristics in the temperature range of -10°C to +45°C.		Fingerprint powder	
	Magnetic Latent Fingerprint Powder 1 container for each type per kit	Types: - black double contrast - fluorescent orange - fluorescent green and fluorescent red Weight: - 200 gr per container		Magnetic Latent Fingerprint Powder 1 container for each type per kit	
	Ostrich/Marabu Brush	Type: - Transparent plastic anti-static sprayer, individually packed. Working length of the brush: - length of at least 80 mm and a thickness at the - thickness at the handle joint of at least 8 mm		Ostrich/Marabu Brush	
	Squirrel/Camel Hair Brushes Minimum Quantity: 3 pieces	Shape: Semi-Round Trim: Untrimmed Length: 14 cm and 19 cm (respectively for the two sizes)		Squirrel/Camel Hair Brushes Minimum Quantity: 3 pieces per set	

	per set	Handle Material: Wooden Handle Length: Standard Active Part Length: 30-40 mm Active Part Thickness: At least 8 mm at the joint with the handle Packaging: Individual packaging for each brush set			
	Glass fiber wire brush: Minimum 4 pieces per set	Active Part Length: Minimum 50 mm Active Part Thickness: Minimum 8 mm at the joint with the handle Handle Material: Plastic or Wooden Handle Length: Standard		Glass fiber wire brush: Minimum 4 pieces per set	
	Black anodised aluminum magnetic brush 2 pcs per kit	Equipped with spring rod for demagnetization.		Black anodised aluminum magnetic brush 2 pcs per kit	
	Master kit small particle reagent (SPR): 2 containers with pump, 250-500 ml	Description: Set of reagents for revealing latent papillary marks on wet, white and black surfaces.		Master kit small particle reagent (SPR): 2 containers with pump, 250-500 ml	
	Gelatine Lifters/ Gelatine sheets for the lift of latent papillary mark: min. 10 sheets of each type per kit	Dimensions: 13x18 cm (±2 cm) Type: white, black and transparent Components: <u>Support:</u> <ul style="list-style-type: none"> ● Material: Plastic or paper laminated on both sides ● Properties: Resistant to mechanical stretching, malleable, unaffected by humidity conditions, opaque (for black foil) and transparent (for transparent foil), allows writing without affecting adhesive 		Gelatine Lifters/ Gelatine sheets for the lift of latent papillary mark: min. 10 sheets of each type per kit	

		<u>Adhesive Component:</u> <ul style="list-style-type: none"> ● Homogeneous adhesive substance, uniformly distributed over the entire surface ● Adhesive to powders used to reveal papillary marks and protective film ● Adhesion to support without absorption or diffusion of powders ● Stable physical-chemical characteristics (-10°C to +45°C), black (for black foil) and transparent (for transparent foil) ● Allows gluing, removal, and regluing of protective film <u>without affecting adhesive</u> <u>Protective Layer:</u> <ul style="list-style-type: none"> ● Material: Transparent plastic foil, malleable, no adhesion to embossing powder ● Easily removable from the adhesive surface. 			
	Instant adhesive sheets for lifting latent papillary mark: min. 10 sheets for each type per ki	Dimensions: 10x15 cm (±2 cm) Type: White, black and transparent.		Instant adhesive sheets for lifting latent papillary mark: min. 10 sheets for each type per ki	
	Latent papillar mark lift foil set. transfer film and fixation of traces:	Dimensions: 5x10 cm (±2 cm) Type: White, black and transparent.		Latent papillar mark lift foil set. transfer film and fixation of traces: min 50 sheets for each type per kit	

	min 50 sheets for each type per kit				
	Adhesive film with mm and cm measuring unit: 2 pcs, each 20 m of each type per kit.	Type: White, black and transparent. Width: 10 - 14 mm		Adhesive film with mm and cm measuring unit: 2 pcs, each 20 m of each type per kit.	
	Case:	Metal cover with tight magnetic closure		Case:	
	Rubber roller:	Width of the active surface: 5-7 cm Handle: Wood		Rubber roller:	
	Fingerprint Ink 1 container per kit Pad:	Solution/gel for cleaning mascara, 1 container of at least 100 ml;		Fingerprint Ink 1 container per kit Pad:	
	Hemident reagent 1 bottle per kit			Hemident reagent 1 bottle per kit	
	Stainless steel scissors			Stainless steel scissors	
	Plastic container for brushes and powder:	Durable plastic material to withstand various environmental conditions and dust		Plastic container for brushes and powder:	
	Flat tweezers:	Material: stainless steel Length: 23-26 cm		Flat tweezers:	
	Sharp tweezers:	Material: stainless steel Length: 8-10 cm		Sharp tweezers:	
	Device for cutting glass:			Device for cutting glass:	
	Magnifying glass :	Diameter: 75 mm Magnification: not less than 3.5x		Magnifying glass :	
	Metal measuring tape:	Length: 3m		Metal measuring tape:	
	Ink cleaning solution: 1 bottle per kit	Volume: 250 ml or bigger		Ink cleaning solution: 1 bottle per kit	
	Flashlight with clamp:	Illumination: - High-intensity light source, providing bright illumination for		Flashlight with clamp:	

	<p>various applications</p> <ul style="list-style-type: none"> - Equipped with efficient LED bulbs <p>Clamp Mechanism:</p> <ul style="list-style-type: none"> - The flashlight incorporates a sturdy clamp mechanism for secure attachment to surfaces or objects. - The clamp is designed to provide firm grip and stability, allowing for hands-free operation in diverse environments. - Adjustable clamp size to accommodate a wide range of objects, such as poles, handles, or surfaces of varying thickness. <p>Power Source:</p> <ul style="list-style-type: none"> - Powered by rechargeable batteries. - Offers convenient charging options, such as USB charging ports or docking stations, for hassle-free recharging. 			
Transparent lifting tape:	Minimum length: 15 m Width: 5-9 cm		Transparent lifting tape:	
Transparent lifting tape:	Minimum length: 15 m Width: 3-5 cm		Transparent lifting tape:	
Ninhydrin Spray: 1 bottle	Volume: 100 ml or more		Ninhydrin Spray: 1 bottle	
Rhodamine Crystals 6G	Container Type: Glass container with threaded cap Capacity: 25 grams		Rhodamine Crystals 6G	
Gentian Violet Crystals:	Container Type: Glass container with threaded lid Capacity: 25 grams		Gentian Violet Crystals:	

	1,8-Diazafluoren-9-one bottle per kit	Type: solution Container Type: Glass or plastic container with screw cap Capacity: 120 milliliters		1,8-Diazafluoren-9-one bottle per kit	
	Portable lamp with UV rays:	Power Output: Minimum equivalent power of 6 watts. Wavelength: Emitting UV light with a wavelength of 365 nanometers (nm), suitable for specific fluorescence excitation. Dimensions: Approximately 7-10 cm in width, 30-35 cm in height, and 5-8 cm in depth, providing a compact and portable form factor for easy handling and transportation. Handle: Equipped with a handle for comfortable grip and convenient carrying during use. Power Supply: Operates on either 220 volts AC power supply or battery power. Battery: Includes a rechargeable battery for portable operation. Charger: equipped.		Portable lamp with UV rays:	
	Documentation:			Documentation:	
		User manuals for the kit (if applicable), preferable in Romanian language (optional English or Russian).			